

# Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/IL05/000303

International filing date: 17 March 2005 (17.03.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US  
Number: 60/553,965  
Filing date: 18 March 2004 (18.03.2004)

Date of receipt at the International Bureau: 09 June 2005 (09.06.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland  
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

PCT/IL 2005 / 000303

25 MAY 2005

PA 1313209

**THE UNITED STATES OF AMERICA**

**TO ALL TO WHOM THESE PRESENTS SHALL COME:**

**UNITED STATES DEPARTMENT OF COMMERCE**

**United States Patent and Trademark Office**


**May 02, 2005**

**THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM  
THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK  
OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT  
APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A  
FILING DATE UNDER 35 USC 111.**

**APPLICATION NUMBER: 60/553,965**

**FILING DATE: *March 18, 2004***

**By Authority of the  
COMMISSIONER OF PATENTS AND TRADEMARKS**



*H. L. Jackson*  
**H. L. JACKSON**

**Certifying Officer**



13281 U.S. PTO

PTO/SB/16 (01-04)

Approved for use through 07/31/2006. OMB 0651-0032

U.S. Patent and Trademark Office: U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

**PROVISIONAL APPLICATION FOR PATENT COVER SHEET**

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No.

U.S. PTO  
60/553965

031804

INVENTOR(S)					
Given Name (first and middle [if any])		Family Name or Surname		Residence (City and either State or Foreign Country)	
ELAN		ZIV MD		RAMAT GAN, ISRAEL	
Additional inventors are being named on the _____ separately numbered sheets attached hereto					
TITLE OF THE INVENTION (500 characters max)					
A DISPOSABLE DEVICE FOR TREATMENT & PREVENTION OF PELVIC ORGAN PROLAPSE					
Direct all correspondence to: CORRESPONDENCE ADDRESS					
<input type="checkbox"/> Customer Number: _____					
OR					
<input type="checkbox"/> Firm or Individual Name		ELAN ZIV MD			
Address		8 HAILANOT ST			
Address					
City		RAMAT GAN		State	
Country		ISRAEL		Zip	52648
		Telephone		Fax	+972-3-7369263
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification Number of Pages		8		<input type="checkbox"/> CD(s), Number _____	
<input checked="" type="checkbox"/> Drawing(s) Number of Sheets		5		<input type="checkbox"/> Other (specify) _____	
<input checked="" type="checkbox"/> Application Data Sheet. See 37 CFR 1.76					
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT					
<input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.				FILING FEE Amount (\$)	
<input checked="" type="checkbox"/> A check or money order is enclosed to cover the filing fees.				<div style="border: 1px solid black; padding: 10px; display: inline-block;">\$80</div>	
<input type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: _____					
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.					
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.					
<input checked="" type="checkbox"/> No.					
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____					

[Page 1 of 2]

Respectfully submitted,  
SIGNATURE [Signature]  
TYPED or PRINTED NAME D.R.E. ZIV MD  
TELEPHONE +972-53-569928

Date 15.3.04  
REGISTRATION NO. \_\_\_\_\_  
(if appropriate)  
Docket Number: \_\_\_\_\_

**USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT**

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

13281  
U.S. PTO

PTO/SB/17 (10-03)

Approved for use through 07/31/2009. OMB 0851-0032  
U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

# FEE TRANSMITTAL for FY 2004

Effective 10/01/2003. Patent fees are subject to annual revision.

☒ Applicant claims small entity status. See 37 CFR 1.27

TOTAL AMOUNT OF PAYMENT (\$ 80

## Complete if Known

Application Number  
Filing Date MARCH 15, 2004  
First Named Inventor Dr Elan Ziv  
Examiner Name  
Art Unit  
Attorney Docket No.

## METHOD OF PAYMENT (check all that apply)

☐ Check ☐ Credit card ☒ Money Order ☐ Other ☐ None

☐ Deposit Account:

Deposit  
Account  
Number  
Deposit  
Account  
Name

The Director is authorized to: (check all that apply)

☐ Charge fee(s) indicated below ☐ Credit any overpayments

☐ Charge any additional fee(s) or any underpayment of fee(s)

☐ Charge fee(s) indicated below, except for the filing fee  
to the above-identified deposit account.

## FEE CALCULATION

### 1. BASIC FILING FEE

Large Entity Fee Code	Small Entity Fee Code	Fee Description	Fee Paid
1001 770	2001 385	Utility filing fee	
1002 340	2002 170	Design filing fee	
1003 530	2003 265	Plant filing fee	
1004 770	2004 385	Reissue filing fee	
1005 160	2005 80	Provisional filing fee	80

SUBTOTAL (1) (\$ 80

### 2. EXTRA CLAIM FEES FOR UTILITY AND REISSUE

Total Claims	Extra Claims	Fee from below	Fee Paid
Independent Claims	-20** =	X	
Multiple Dependent	-3** =	X	

Large Entity Fee Code	Small Entity Fee Code	Fee Description
1202 18	2202 9	Claims in excess of 20
1201 86	2201 43	Independent claims in excess of 3
1203 280	2203 145	Multiple dependent claim, if not paid
1204 86	2204 43	** Reissue independent claims over original patent
1205 18	2205 9	** Reissue claims in excess of 20 and over original patent

SUBTOTAL (2) (\$

\*\*or number previously paid, if greater; For Reissues, see above

## FEE CALCULATION (continued)

### 3. ADDITIONAL FEES

Large Entity Small Entity

Fee Code	Fee (\$)	Fee Code	Fee (\$)	Fee Description	Fee Paid
1051 130		2051 65		Surcharge - late filing fee or oath	
1052 50		2052 25		Surcharge - late provisional filing fee or cover sheet	
1053 130		1053 130		Non-English specification	
1812 2,520		1812 2,520		For filing a request for ex parte reexamination	
1804 920*		1804 920*		Requesting publication of SIR prior to Examiner action	
1805 1,840*		1805 1,840*		Requesting publication of SIR after Examiner action	
1251 110		2251 55		Extension for reply within first month	
1252 420		2252 210		Extension for reply within second month	
1253 950		2253 475		Extension for reply within third month	
1254 1,480		2254 740		Extension for reply within fourth month	
1255 2,010		2255 1,005		Extension for reply within fifth month	
1401 330		2401 165		Notice of Appeal	
1402 330		2402 165		Filing a brief in support of an appeal	
1403 280		2403 145		Request for oral hearing	
1451 1,510		1451 1,510		Petition to institute a public use proceeding	
1452 110		2452 55		Petition to revive - unavoidable	
1453 1,330		2453 665		Petition to revive - unintentional	
1501 1,330		2501 665		Utility issue fee (or reissue)	
1502 480		2502 240		Design issue fee	
1503 640		2503 320		Plant issue fee	
1460 130		1460 130		Petitions to the Commissioner	
1807 50		1807 50		Processing fee under 37 CFR 1.17(q)	
1808 180		1808 180		Submission of Information Disclosure Stmt	
8021 40		8021 40		Recording each patent assignment per property (times number of properties)	
1809 770		2809 385		Filing a submission after final rejection (37 CFR 1.129(a))	
1810 770		2810 385		For each additional invention to be examined (37 CFR 1.129(b))	
1801 770		2801 385		Request for Continued Examination (RCE)	
1802 900		1802 900		Request for expedited examination of a design application	

Other fee (specify)

\*Reduced by Basic Filing Fee Paid

SUBTOTAL (3) (\$

## SUBMITTED BY

Name (Print/Type) Dr Elan Ziv MD

Signature

Registration No.  
(Attorney/Agent)

(Complete if applicable)

Telephone +97237369261

Date March 14, 2004

**WARNING: Information on this form may become public. Credit card information should not be included on this form. Provide credit card information and authorization on PTO-2038.**

This collection of information is required by 37 CFR 1.17 and 1.27. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 12 minutes to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

## **A DISPOSABLE DEVICE FOR TREATMENT & PREVENTION OF PELVIC ORGAN PROLAPSE.**

The present invention relates generally to the field of treatment and prevention of pelvic organ prolapse in female patients. The invention describes a vaginal disposable device which is inserted and removed in a no-self-touch technique, by the patient herself, using a disposable applicator.

**Inventor:** Dr Elan Ziv, MD OBGYN, Urogynecologist

### **Background of the Invention**

#### **General Data**

##### **Definition**

Pelvic organ prolapse is defined as a condition in which vaginal wall support is lost, and various pelvic organs prolapse into the vagina. This is a very bothering condition, though in most cases it is not a dangerous one. It might appear alone or in combination with urinary stress incontinence.

POP (Pelvic Organ Prolapse) is a very common condition in females, in which various organs, surrounding the vagina, prolapse into it. The reasons for such prolapse are numerous, mainly because the vagina has a very low tensile strength within its walls due to damage to muscles, nerves, fascias etc. There might also be a change within the collagen content, thereby causing a weaker pelvic floor. According to the older definition, prolapse may be divided into five categories, according to the organ that is sagging down (urethra, bladder, uterus, rectum and the pouch of Douglas (small bowl)), and to three grades according to the amount of descent (within the vagina, at the entrance to the vagina, protruding out of the vagina). There might be a combination of various organs prolapse at the same time, with different level of descent. The newer classification (POP-Q) takes into account other factors, such as location of the prolapse and the distance from the entrance of the vagina. There is very little data regarding the prevalence of the problem but it seems to be age related, and older women have much higher tendency towards developing some form of pelvic organ prolapse.

### High risk group for Pelvic Organ Prolapse

There are very few statistical & Epidemiological studies regarding the causes for the development of pelvic organ prolapse and high risk groups. The main risk factors are:

- Advanced age
- Pregnancy
- Delivery. During labor, there are many tears, stretchings and hematomas within the muscles, the ligaments and the nerves of the pelvic floor. Instrumental delivery might be a further factor.
- Number of deliveries
- Race. White women tend to have more prolapses than black or oriental women.
- Menopause
- Obesity
- Hard labor.

### Prevalence

It was found that the prevalence of severe pelvic organ prolapse in women with jobs in labor/factory positions was 11.1%, 6.6% among those working at home, 3.0% in professional/ managerial jobs, 3.3% in service, and 0.8% among technical/sales/clerical workers. Laborers/factory workers had significantly more severe prolapse than the other job categories ( $P < .001$ ). Table 1 summarizes some article on the subject

APPENDIX I – Prevalence Rates of POP				
	Samuelsson et al 1999	Bland et al 1999	Progetto Menopausa 2000	Swift 2000
Country	Sweden	U.S.	Italy	U.S.
N	487	241	21,449	497
Age range	20-59	45-55	Around menopause I	8 - >70
Race				
White	100%	89%	100%	47%
Black		9%		52%
Parous	54%	Unspecified	85.6%	93%
Hysterectomy	4%	28%	None	28%
Staging technique	Above, to, beyond hymen	ICS	Uterus only Baden 0, 1, ≥2	ICS
Stage				
0	71.2%	73%	94.5%	6.4%
1	28.8%	23%	3.6%	43.3%
2	2%	4%	1.9%	47.7%
3	0	0		2.6%
4	0	0		0
Source of subjects	76% of every 3 year screen participants	Responded to ad for soybean supplement study	268 first-level menopause clinics in Italy.	Routine gynecologic examination

TABLE 1

"A disposable device for treatment & prevention of pelvic organ prolapse "  
Dr Elan Ziv, MD OBGYN, Urogynecologist

### **Anatomy and pathophysiology**

Understanding the female genital anatomy is extremely important in order to understand the basis of this invention.

FIG. 1A is an anatomical frontal cut section of the vagina

5 FIG. 1B is a front view of vagina creating the shape of the letter H.

FIG. 1C is a side view of the vaginal diameters.

FIG. 2 is a side view of the female pelvis.

The vagina is a hollow organ with multiple shapes and diameters at different depths.

10 The sidewalls (FIG.1A) have thin longitudinal and circular muscles (8), externally wrapped by the pelvic fascia (10), which tends to become weakened, stretched and attenuated over the years and allow for sagging and prolapse of different pelvic organs into and through the vagina. Internally the vagina is covered with the mucosa (6).

Although potentially being a tube, the vaginal cavity is almost always obliterated with  
15 the anterior wall (12) and the heavy bladder on it collapse and fall into the vaginal lumen. The posterior wall is pushed anteriorly by the bowel (14), thereby creating the shape of the letter H (FIG.1B). Vaginal diameters, either lateral or antero-posterior (FIG.1C, FIG.2) are never the same along vaginal axis. Lateral diameter is usually shorter closer to the entrance (18) gradually becoming longer internally (16). The same  
20 applies for the antero-posterior diameter, while diameter close to the entrance (44) is usually shorter than the one behind the perineal body (48). The vagina may therefore be looked at as a funneled shape organ.

There are certain defined organs in direct contact with the vagina or within short proximity. The bladder (36) is a hollow, sack like organ, in which the urine is  
25 accumulated prior to it's expulsion outside the body through the urethra (42). The bladder is located behind the pubic bone (38) and rests on the middle third of the vagina. The urethra is a short muscular pipe, 25-35 mms in length, in direct contact with the bladder at the bladder neck (40) resting on the lower third of the anterior vaginal wall, in direct contact with the pubic bone. The uterus (30) is a pear shaped organ that  
30 has a lower part, the cervix (32) which protrudes into the dome of the vagina. The body of uterus is an abdominal organ, which, in most cases bands forward on the bladder. The cervix, while protruding into the vagina creates 4 fornices – 2 laterals which are of the same dimensions, one smaller anterior (34), and one posterior (22), which is actually

behind the cervix. The posterior fornix has direct contact with the pouch of Douglas (20) which is the deepest part of the abdominal cavity, on top of the vagina. The posterior surface of the vagina is in contact with the rectum (24), until the area where a thick muscular body, the perineal body (26) which is part of the pelvic floor and the sphincter of the anus (28).

Prolapse of adjacent organs into the vagina is quite common, in variable degrees. The reasons for such a prolapse are mainly damage to the endopelvic fascia which surrounds the organs and keeps them in the right position, pelvic floor muscular damage, and neural damage. Prolapse of the urethra, named urethrocel, is one of the several possibilities for pelvic organ prolapse around the vagina. The most frequent prolapse is of the bladder, named cystocele, as a result of its weight. Prolapse of posterior fornix, with bowel inside is called enterocele, and prolapse of the posterior vaginal wall with the rectum in it is termed rectocele. When the ligaments that hold the uterus in place weaken, a uterine descent occurs. In cases where the uterus has been removed, and the vagina is dome shaped, a vault prolapse may occur.

### **Current treatment**

At present – there are 2 ways of dealing with Pelvic Organ Prolapses:

1. Surgical – vaginal or abdominal
2. Use of vaginal devices (pessaries) that are inserted into the vagina and mechanically reduce the prolapse by pushing the wall aside and upwards.

Efforts to avoid surgical procedures have resulted in the development of a number of non-surgical vaginal devices, inserted into the vagina by the surgeon or the patient.

Vaginal devices are well known for their tremendous diversity in shapes and sizes. These devices are meant to prevent prolapse of the vaginal walls, with different location of pressure application.

Some of these devices tend to block all flow of urine from the bladder. Therefore, when a patient needs to urinate, the device must be removed from the vagina or must be collapsed to remove the pressure applied against the bladder neck. Trying to solve this problem, vaginal devices were developed in special shapes, without completely blocking the bladder neck so that the patient may urinate with the device in place ( for example, US patents 6189535 B1, 6158435, 5894842, 5611768). These devices,

however, are generally large and intrusive and, therefore, are uncomfortable to wear, with low patient satisfaction and compliance. They are also relatively expensive, and therefore designed to be reusable.

5 There are several drawbacks of existing vaginal devices:

1. Most of them are intended to be reusable, hence they are made as resilient large bodies, made of plastic, hard rubber, or other such materials, in order to preserve their shape and function for a long time.
- 10 2. Insertion of large noncompliant bodies is sometimes difficult, painful or unpleasant, sometimes necessitating a medical practitioner. Removal causes same unpleasantness.
3. Most reusable devices are meant to remain in the vagina for prolonged periods of time, thereby causing irritations, pressure ulcers, infections, foul smelling discharge, etc.
- 15 4. Such reusable devices have a long standing bad reputation among patients and medical practitioners for being unpleasant, causing infectious discharge and foul odor, and being associated with disability and old age.
- 20 5. Some reusable devices are meant to be inserted daily by the patient, and to be removed after several hours, by means of pulling a string or with a finger, to be cleaned for re-use, and to be kept in certain conditions prior to following insertion. Some patients are reluctant to touch their selves in such intimate parts of their bodies, or disgusted to clean the device, hence their reluctance to use it.

Vaginal devices may cause three main side effects:

- 25 1. **Infections**-any foreign material, anywhere in the body, may become infected by several kinds of organisms, and cause formation of a foul smelly discharge. Reusable devices are certainly prone to cause such infections, but also, less frequently, disposable ones. In order to prevent such infections, devices should be:
  - 30 • disposable
  - used for a limited length of time, in order to prevent organism from growing, and to prevent damage, by vascular pressure, to the vaginal wall.

- made in a way that will allow vaginal or cervical normal secretions to flow out of the body
- made of certain known materials that would not permit growth of organisms.

- 5        2. **Toxic Shock Syndrome (TSS)** is a condition, described some years ago, in which abundant growth of *Staphylococcus Aureus*, a bacteria that utilizes the cellulose within sanitary menstrual tampons, released large amounts of toxin. That toxin caused a collapse of vital body systems, forming a dangerous condition. In order to overcome this, sanitary tampons do not contain cellulose
- 10        anymore, and vaginal devices, as described earlier, should be made of materials that would not permit the growth of bacteria, and allow for discharge flow. This is best done by using properly made disposable devices.
- 15        3. **Pressure Necrosis** – prolonged pressure on the vaginal sidewalls may cause pressure on blood vessels with resultant necrosis, bleeding and infections. This might be prevented by using disposable devices only, for a predetermined length of time.

20

### **The invention**

FIG. 2 is a side view of the female pelvis.

FIG. 3A is a top view of the invention.

FIG. 3B is a side view of the invention.

5 FIG. 4 is a perspective view of the applicator.

FIG.5A is a sectional side view of the device within the applicator.

FIG.5B is a sectional top view of the device within the applicator.

10 The device (FIG. 3) is intended to be applied into the vagina with an applicator much like the insertion of a regular menstrual tampon .After insertion, the device should expand significantly to the predefined shape and size, thereby exerting predefined appropriate pressure on both lateral vaginal walls, pushing them aside. The apex of the vagina shall be pushed upwards at the same time. That expansion of the device shall eventually create linear stretching of the anterior & posterior vaginal walls, while  
15 creating a new shape of intra-vaginal hollow (rectangle). The device may be left inside the vagina for several hours. Removal for disposal will occur while pulling a string and collapsing the device to a much smaller size.

The device (FIG. 3A), from a top view, looks like a round ring (50+52) and a string (54). It is made of elastic materials, such as silicon, polyurethane, etc. This flexible  
20 material gives the device the ability to be folded into the applicator, and to be pushed out easily in small dimensions into the vagina. The device is covered by a sponge rubber (52) which makes it more comfortable to wear. The sponge is also reducing the pressure on the vagina's sidewalls and helps preventing pressure necrosis. FIG.3B is a side view of the ring what shows its convex shape. This unique shape helps the ring to  
25 be folded and inserted easily into the applicator.

As with an ordinary menstrual tampon, (FIG.4) the device will be kept and supplied within an applicator, ready for use. FIG.5A is a sectional side view of the device (56) inside the applicator (58). FIG.5B is a sectional top view of the device (60) inside the applicator (62).

30 Removal of the device is inherited within the device itself. Pulling the string (54), which will protrude outside of the vagina as with the regular menstrual tampon, will cause the device to reduce in size and be pulled out of the vagina by the same pull of the string, in order to be wasted.

Discharge will be allowed to flow out of the vagina, as due to its shape, the device does not block the vagina.




5 The device is intended for use in women who suffer from a prolapse of at least one of the above mentioned organs. Preferably, it should aim at women who suffer from any prolapse descending as low as the level of the hymen, but it might just as well be used in women who have greater degrees of prolapse

The invention has in its basic concept the following features:

- Being a disposable device.
- 10 • Insertion of the device is always with an applicator.
- Easy & comfortable insertion and removal.
- Being comfortable to wear.
- Being hygiene & odorless
- 15 • Being a familiar procedure to most female patients – as inserting a menstrual tampon.
- Being inserted by the patient herself, in a no-self-touch technique, with a disposable inserter.
- Being removed by the patient herself, in a no-self-touch technique, with the device collapsing and becoming of small size for painless removal.
- 20 • Being of high availability, easy to get everywhere, sold as an Over the Counter (OTC) device.
- Being of low cost.
- Having complete confidentiality, as with the use of menstrual tampons.
- Having the ability to be removed instantly when needed.
- 25 • No blockage of vaginal discharge.

**Alternative embodiments of the invention.**

- It may be manufactured in different sizes
- 30 • It may be made of many flexible materials, such as plastics, cardboard, etc.








Fig.2


"A disposable device for treatment & prevention of pelvic organ prolapse"  
 Dr Elan Ziv, MD OBGYN, Urogynecologist

BEST AVAILABLE COPY



"A disposable device for treatment & prevention of pelvic organ prolapse"  
 Dr Elan Ziv, MD OBGYN, Urogynecologist

BEST AVAILABLE COPY



"A disposable device for treatment & prevention of pelvic organ prolapse "  
Dr Elan Ziv, MD OBGYN, Urogynecologist

BEST AVAILABLE COPY




Fig. 5A




Fig. 5B

"A disposable device for treatment & prevention of pelvic organ prolapse"  
 Dr Elan Ziv, MD OBGYN, Urogynecologist

BEST AVAILABLE COPY